



Dropper

Cardiac Markers Control

- Easy-to-use, dropper bottles
- Liquid, made with real human serum
- 90 day unopened and 60 day open vial stability when stored at 2-8 C*
- 3 year frozen stability; can be refrozen 2 times*

Analytes*

CK-MB
CK, Total

CRP
Digoxin

L-Homocysteine
Myoglobin

Troponin I
Troponin T

**Refer to the package insert of currently available lots for specific analyte and stability claims*

Take convenience to a whole new level. We craft Dropper Cardiac Markers Control to monitor the performance of serum cardiac marker assay methods by dropping the control directly into the instrument sample cup with an easy-to-use, dropper bottle for outrageously simple dispensing and maximum use. Dropper Cardiac Markers Control is designed for use with most major chemistry analyzers including Advia®, Access®, Dimension®, Elecsys®, Olympus®, Stratus®, Synchron®, and Vitros®. No hassle, no waste, no problem.

ORDERING INFORMATION:

1011-61	Level 1.....	6 x 3 mL
1012-61	Level 2.....	6 x 3 mL
1013-61	Level 3.....	6 x 3 mL
1010-02	3 Level Set.....	6 x 3 mL



For More Information: phone +1.310.536.0006 email info@4qc.com web www.4qc.com