



# *Dropper Plus*

## Point-of-Care Urine Dipstick Control

- 1 month open vial stability at room temperature for maximum portability\*
- Liquid, made with real human urine
- 18 month unopened and open vial stability when stored at 2-8 C\*

### Analytes\*

Bilirubin	Creatinine	hCG	Leukocytes	Nitrite	Protein	Urobilinogen
Blood	Glucose	Ketones	Microalbumin	pH	Specific Gravity	

*\*Refer to the package insert of currently available lots for specific analyte and stability claims*

Take your control with you. We craft Dropper Plus POC Urine Dipstick Control to make point-of-care urine dipstick quality control simple. One month of room temperature stability eliminates storage problems and provides maximum portability – perfect for sites without refrigeration. And easy-to-use, dropper bottles make dispensing outrageously simple. We include values for all dipstick analytes plus microalbumin and creatinine; and qualitative results for hCG early pregnancy detection test methods. Dropper Plus POC Urine Dipstick Control is designed for use with most reagent strips. It can also be used for confirmatory tests and refractometry. Have control - will travel.

### ORDERING INFORMATION:

1440-04 2 Level Set.....10 x 5 mL



**For More Information:** phone +1.310.536.0006 email [info@4qc.com](mailto:info@4qc.com) web [www.4qc.com](http://www.4qc.com)