



# Dropper

## Spinal Fluid Control

- Easy-to-use, dropper bottles
- Liquid, human protein matrix
- 2 year unopened and 90 day open vial stability when stored at 2-8 C\*

### Analytes\*

Albumin	Glucose	Lactic Acid	Sodium
Chloride	IgG	LDH	Total Protein

*\*Refer to the package insert of currently available lots for specific analyte and stability claims*

Take convenience to a whole new level. We craft Dropper Spinal Fluid Control to monitor the performance of spinal fluid chemistry assay methods by dropping the control directly into the instrument sample cup with an easy-to-use, dropper bottle for outrageously simple dispensing and maximum use. Dropper Spinal Fluid Control is designed for use with most major chemistry analyzers and nephelometers: Synchron®, Dimension®, AU® series, Vitros®, Hitachi®, Integra®, and Pyrogallol Red microprotein methods. No hassle, no waste, no problem.

### ORDERING INFORMATION:

- 1451-31 Level 1.....3 x 3 mL
- 1452-31 Level 2.....3 x 3 mL



**For More Information:** phone +1.310.536.0006 email [info@4qc.com](mailto:info@4qc.com) web [www.4qc.com](http://www.4qc.com)