



# Dropper

## Urine Dipstick Control

- Easy-to-use, dropper bottles
- Liquid, made with real human urine
- 18 month unopened and open vial stability when stored at 2-8 C\*

### Analytes\*

Bilirubin	Creatinine	hCG	Leukocytes	Nitrite	Protein	Urobilinogen
Blood	Glucose	Ketones	Microalbumin	pH	Specific Gravity	

*\*Refer to the package insert of currently available lots for specific analyte and stability claims*

Take convenience to a whole new level. We craft Dropper Urine Dipstick Control to monitor the performance of visual and instrumental readings of urine dipsticks by dropping the control directly onto the dipstick with an easy-to-use, dropper bottle for outrageously simple dispensing and maximum use. We include values for all dipstick analytes plus microalbumin and creatinine; and qualitative results for hCG early pregnancy detection test methods.

Dropper Urine Dipstick Control is designed for use with most reagent strips. It can also be used for confirmatory tests and refractometry. No hassle, no waste, no problem.



### ORDERING INFORMATION:

1440-02 2 Level Set.....4 x 25 mL

**For More Information:** phone +1.310.536.0006 email [info@4qc.com](mailto:info@4qc.com) web [www.4qc.com](http://www.4qc.com)