



Microbumin

Microalbumin Control

- 3 year unopened and 6 month open vial stability when stored at 2-8 C*
- Liquid, made with real human urine
- Includes creatinine

Analytes*

Albumin

Creatinine

**Refer to the package insert of currently available lots for specific analyte and stability claims*

Extraordinary stability makes it super convenient. We craft Microbumin Microalbumin Control to monitor the performance of popular strips including Clinitek® Microalbumin Reagent Strips and Micral™ test strips. Microbumin is also perfect for use with turbidimetric immunoassays using Advia®, Clinitek®, DCA 2000®, Dimension®, Integra®, Olympus®, Synchron®; and nephelometry. Longer stability, less waste.



ORDERING INFORMATION:

1250-01 2 Level Set.....10 x 7 mL

For More Information: **phone** +1.310.536.0006 **email** info@4qc.com **web** www.4qc.com