



QuanTscopics®

Urine Microscopics Control

- Cost effective
- Liquid, made with real human urine
- 2 year unopened and open vial stability when stored at 2-8 C*

Analytes*

Crystals

Red Blood Cells

White Blood Cells

**Refer to the package insert of currently available lots for specific analyte and stability claims*

Focus on the cells. We craft QuanTscopics to give you a simple, cost effective way to monitor the performance of urine microscopic analysis. And easy-to-use, squeeze bottles make dispensing outrageously simple. We include qualitative and quantitative values for standardized urinalysis systems and analyzers. Take a closer look and see what you've been missing.



ORDERING INFORMATION:

1481-41 Level 1..... 4 x 120 mL
 1482-41 Level 2..... 4 x 120 mL

For More Information: phone +1.310.536.0006 email info@4qc.com web www.4qc.com