



Spinalscopics™

Spinal Fluid Cell Count Control

- 1 month open vial stability at room temperature*
- Liquid, human protein matrix with stabilized human cells
- 18 month unopened and 6 month open vial stability when stored at 2-8 C*

Analytes*

Red Blood Cells

White Blood Cells

**Refer to the package insert of currently available lots for specific analyte and stability claims*

Focus on the cells. We craft Spinalscopics Spinal Fluid Cell Count Control to monitor the performance of red and white blood cell count methods in cerebrospinal fluid samples performed manually using a hemocytometer. And we provide an easy-to-use, dropper dispenser for outrageously simple use. Take a closer look and see what you've been missing.

ORDERING INFORMATION:

- 1461-31 Level 1.....3 x 3 mL
1462-31 Level 2.....3 x 3 mL



For More Information: [phone +1.310.536.0006](tel:+13105360006) [email info@4qc.com](mailto:info@4qc.com) [web www.4qc.com](http://www.4qc.com)